

# Comfortmax Elite Cervical Orthosis

Designed for immobilization of the cervical spine, arthritis, spinal stenosis, and post trauma use, the Comfortmax Elite Cervical Orthosis features an aluminum, multi-post anterior adjustment system, allowing full linear adjustment. The CV-500 is lightweight and vented; promoting air flow, improving patient comfort, and making it MRI compatible. The unique strapping design allows the Elite to fit a wide variety of patients, eliminating the need for multiple sizes.




## Universal Sizing

The CV-500's unique design fits most patients easily, eliminating the need to stock multiple sizes



## Easy Application and Removal

A hook-and-loop system allows for easy donning/doffing without requiring assistance



## Total Immobilization

The CV-500 features a sturdy, aluminum anterior yoke, providing optimum support and immobilization to help prevent further injury